
Life Sciences

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	Acetic acid, glacial	Acetic acid, 90%	Acetic acid, 30%	Acetic acid, 10%	Hydrochloric acid, conc. (35%)	Hydrochloric acid, 6N, (20%)	Nitric acid, conc. (3.3%)	Nitric acid, 6N (67%)	Sulfuric acid, conc. (96%)	Sulfuric acid, 6N (16%)	Amyl alcohol	Benzyl alcohol	Butanol	Ethanol	Isopropanol	Methanol	Ammonium hydroxide, 3N (5.7%)	Ammonium hydroxide, 6N (11.4%)	Potassium hydroxide, 3N (15%)	Sodium hydroxide, 3N (11%)	Amyl acetate	Butyl acetate	Cellosolve acetate	Ethylacetate	Isopropyl acetate	Methyl acetate	Ethyl ether	Tetrahydrofuran	Ethylene glycol/water(50/50,v/v)	Glycerol	Propylene glycol	Benzene	Xylene	Carbon tetrachloride	Chloroform	Ethylene dichloride	Methylene chloride	Tetrachloroethylene	Acetone	Cyclohexanone	Methyl ethyl ketone (MEK)	Methyl isobutyl ketone	Cottonseed	Peanut	Acetonitrile	Dimethyl formamide (DMF)	Dimethyl sulfoxide (DMSO)	Formaldehyde, 34%	Formaldehyde, 4%	Hexane, dry	Kerosene	Pyridine	18 Megohm water
	Acids		Alcohols		Bases		Esters		Ethers		Glycols		Aromatic Hydrocarbons		Halogenated Hydrocarbons		Ketones		Oils		Miscellaneous																																
GH Polypore (GHP) (PP) Membrane	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R																
Teflo, TF (PTFE) Membrane	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R																
FP Vericel™ (PVDF) Membrane	R	R	R	R	R	R	R	R	N	-	R	R	R	R	R	R	N	N	N	N	R	R	L	R	R	R	R	N	-	R	R	R	R	R	R	R	R																
Nylaflow™ (Nylon) Membrane	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N	N	R	N	R	R	R	R	R	R	R	R	R	R	L	R	R	R	R	R	R	R	R																
Supor® (PES) Membrane	R	R	R	R	R	R	R	N	-	N	-	N	N	R	R	R	R	R	R	R	R	R	R	R	R	R	R	N	-	R	R	R	R	R	R	R	R																

Caution

Alcohol residues that are allowed to dry on a filter may cause stress cracks. Pall Corporation recommends that filters used in alcohol processing should remain alcohol wet, or should be flushed with copious quantities of water to remove residuals prior to drying and subsequent reuse.

R = Resistant

No significant change was observed in flow rate or bubble point of the membrane, was detected no visible indication of chemical attack.

L = Limited Resistance

Moderate changes in physical properties or dimensions of the membrane were observed. The filter may be suitable for short term, non-critical use. Hardware or housing may be suitable for short-term exposure at low pressures and ambient temperatures.

N = Not Resistant

The membrane or housing is basically unstable and is not recommended for use.

- = Insufficient Data

Information is not available. Trial testing is recommended.

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
Membrane	Proteinaceous	General Aqueous	Non-Aggressive Organic	Aggressive Organic	Hydrophobic vs Hydrophilic
GHP	++	++	++	++	Hydrophilic
PTFE	-	-	++	++	Hydrophobic
PVDF	++	++	+	-	Hydrophilic
Nylon	+*	++	+	-	Hydrophilic

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BIO-CHECK LABORATORIES LTD.